




# Introduction to ChIP Sequencing and QC

---

Sharon Freshour, PhD  
St. Jude Children's Research Hospital  
December 11, 2024

# The DNB Bioinformatics Core Team

---



**Cody Ramirez, PhD**

Senior Bioinformatics Research Scientist  
Core Director  
Boston, Massachusetts



**Antonia Chroni, PhD**

Senior Bioinformatics Research Scientist  
New York, New York



**Asha Jacob Jannu, PhD**

Bioinformatics Research Scientist  
Indianapolis, Indiana



**Sharon Freshour, PhD**

Bioinformatics Research Scientist  
St. Louis, Missouri




# Introduction to ChIP Sequencing and QC Workshop Overview

- ChIP-Seq overview
  - What is ChIP-Seq?
  - How is ChIP-Seq done?
  - What are other methods for profiling regulatory elements?
- Submitting samples for sequencing and analysis
- CAB's ChIP-Seq pipelines
  - AutoMapper, QC, and Peak Calling
- Example of ChIP-Seq QC and Peak Calling report from CAB



# ChIP Sequencing Overview

# ChIP sequencing is a useful tool for understanding complex transcriptional regulation



Adapted from [The ENCODE Project Consortium \(2011\). PLOS Biology](#)



# How does ChIP-Seq work?

---

- Uses a combination of chromatin immunoprecipitation (ChIP) and NGS (seq)
  - Antibody selection for proteins of interest
  - Next generation sequencing
  - Assays protein-DNA binding in vivo, across genome
- Complements gene expression profiling, DNA accessibility methods
- Caveats
  - Qualitative, not quantitative, profiles enrichment
  - Need quite a bit of material for standard ChIP-Seq
  - Heterogeneity can be hard to capture
  - Must have good antibodies for selection step





# ChIP-Seq library preparation considerations

## Must have sufficient starting material

- At minimum,  $10^7$  cultured cells recommended for single ChIP experiment


## Success dependent on antibody selection

- Should be specific
- Can monoclonal or polyclonal


## Two general ChIP procedures to choose from:

- Native ChIP (N-ChIP)
- Cross-linking ChIP (X-ChIP)

Polyclonal antibodies



Monoclonal antibodies




Credit: [www.abcam.com/chip](http://www.abcam.com/chip)




# Controls are necessary for ChIP-Seq experiments

- Noise in ChIP experiments not uniform
  - Affected by chromatin confirmation, local biases, mappability
- Need to model background noise to distinguish true peaks
- Input controls are necessary to estimate noise
  - Cross-linked, fragmented DNA without antibody enrichment
  - Recommended one for every immunoprecipitation done
  - If constraints, one per sample group can be sufficient
- Isotype (IgG) controls can also be used
  - Immunoprecipitation with an isotype-matched control
  - Similar to experimental antibody, but non-specific binding
- Can also use positive and negative controls, qPCR to check success
  - Positive: check signal in with expected protein binding
  - Negative: check lack of signal in non-enriched region






# Cross-linking proteins and DNA is (often) the first step for ChIP-Seq library preparation



Adapted from [doi.org/10.1038/nrg2641](https://doi.org/10.1038/nrg2641)




# Next step for ChIP-Seq preparation is fragmenting DNA



Adapted from [doi.org/10.1038/nrg2641](https://doi.org/10.1038/nrg2641)




# Fragmented, protein-bound DNA is immunoprecipitated using specific antibodies



Adapted from [doi.org/10.1038/nrg2641](https://doi.org/10.1038/nrg2641)



# Cross-linking is reversed and DNA is purified for sequencing



Adapted from [doi.org/10.1038/nrg2641](https://doi.org/10.1038/nrg2641)




# ChIP-PCR vs. ChIP-chip vs. ChIP-Seq

ChIP-PCR	ChIP-chip	ChIP-Seq
<ul style="list-style-type: none"><li>• Targeted regions of genome</li><li>• Regions known beforehand</li><li>• Cheaper</li><li>• More time efficient</li><li>• qPCR can allow quantitative comparisons</li><li>• qPCR can confirm ChIP successful</li></ul>	<ul style="list-style-type: none"><li>• Whole genome (but can profile specific regions)</li><li>• Microarray-based</li><li>• 30-100 bp resolution typically</li><li>• Requires ~a few micrograms DNA</li><li>• Useful for broad binding</li></ul>	<ul style="list-style-type: none"><li>• Whole genome</li><li>• Next-generation sequencing</li><li>• Single nucleotide resolution</li><li>• Only requires ~10 – 50 ng of DNA</li><li>• Becoming more cost effective</li><li>• Useful for sharp binding</li></ul>



# Alternative methods to ChIP-Seq



Adapted from [doi.org/10.1038/s41596-020-0373-x](https://doi.org/10.1038/s41596-020-0373-x), [activemotif.com/blog-atac-seq](https://activemotif.com/blog-atac-seq)



# Submitting Samples for Sequencing and Analysis

# Submitting ChIP-Seq samples for Hartwell sequencing

- Submit request for Hartwell sequencing via SRM
- Fill out template spreadsheet to submit (or manually file out table in SRM)
  - Sample info for each well of 96 well plate
- Hartwell pipeline is optimized for 10 ng input DNA per sample, 1 ng input accommodated, 100 pg attempted (submit in 52 uL volume)
- Receive email with SRM order confirmation with QR for sample label
- After sequencing, receive email with fastq file paths, basic sequencing QC

	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P
1	Please use this excel to upload samples for <b>Genome Sequencing Service</b> in SRM2 System															
3	Well Location	Sample Name	Is this a Human Derived Sample?	SJ Tissue Bank #	SJUID	Alternative # (lab, cooperative group, etc.)	Submission Material	Xenograft?	Application	Illumina Sequencer	Run Type	Read Length	Molecules Sequenced	Reference Genome	Please specify Reference Genome	User Comments
4	A01								ChIP-seq	NovaSeq	Single End	50 bps	Default for Selected Application			
5	B01						Nucleic Acid for Prep and Sequencing User Made Libraries to be Sequenced									
6	C01															
7	D01															
8	E01															
9	F01															

# Requesting ChIP-Seq analysis from Center for Applied Bioinformatics (CAB)

- CAB AutoMapper for ChIP-Seq runs automatically after Hartwell sequencing completed
- Receive email when AutoMapper submission is started
- Receive email when AutoMapper run is completed
- AutoMapper pipeline does not include analysis, only alignment
- Must submit new SRM request for analysis
- QC and Peak Calling is standard analysis
- Additional analysis considered customized

Select the data type\*

ChIPseq/Cut-and-Run ▾

Select a bioinformatics analysis from the list below. If it's not on our list, please choose "Other" and describe the analysis in the box to the side.


Bioinformatics analysis\*

☒ A1a. QC and Peak calling done  
☐ A1b. Peak annotation (ChIPseq)  
☐ A2. Differential binding analysis (ChIP)  
☐ A3. Chromatin state assignment (chromHMM)  
☐ A14. Competitive Mapping (for spike-in or de-contamination)  
☐ A15a. GO enrichment analysis (ChIP WGBS)  
☐ A16a. DNA motif analysis (ChIPseq/WGBS)  
☐ A17. Super-enhancer analysis and Circular regulatory circuitry (H3K27ac)  
☐ A18. Co-localization analysis  
☐ A26. Spike-in or SpikeInFree normalization  
☐ U2. Mapping (ChIPseq)  
☐ Other




# **CAB ChIP-Seq AutoMapper Pipeline, QC, and Peak Calling**

# CAB AutoMapper pipeline






# CAB QC and Peak Calling

Automapper






# Cross-correlation analysis indicates quality of sequencing, fragment length, binding sites



Adapted from [doi.org/10.1101/gr.136184.111](https://doi.org/10.1101/gr.136184.111)



# Interpreting cross-correlation plots, relative stranded correlation (RSC) values




Qtag	-2	-1	0	1	2
RSC	0, 0.25	0.25, 0.5	0.5, 1	1, 1.5	≥1.5

$$RSC = \frac{cc(fragment\ length) - min(cc)}{cc(read\ length) - min(cc)}$$

Adapted from [doi.org/10.1101/gr.136184.111](https://doi.org/10.1101/gr.136184.111)

# Cross-correlation analysis is only relevant for single-end sequencing

- If sequencing was paired end, ignore cross-correlation results
- Results are meaningless
- Paired reads will have mate read on opposite strand
- Separated by distance of ~average fragment length



# Common quality control metrics for ChIP-seq

- Mapping rate > 80%
- Duplication rate <= 30%
- Total number of unique reads (per ENCODE):
  - > 10M for narrow peaks (point-source data)
  - > 20M for broad peaks
- Fragment size > 100bp
- Qtag > 0 (single-end data only)
  - Wave pattern in cross-correlation plot, RSC invalid
- Visualize data with IGV
  - Check clear peaks
  - Known markers


Qtag	-2	-1	0	1	2
RSC	0, 0.25	0.25, 0.5	0.5, 1	1, 1.5	≥1.5




# CAB pipeline can call broad or narrow peaks

- Peak calling identifies regions with enriched protein-DNA interactions
- Narrow peak calling is done with **MACS2**
- Broad peak calling is done with **SICER**
- If known, can tell CAB broad or narrow
- Or choose auto and pipeline will choose method
- Important to use correct method for peak calling

A



B



Adapted from [doi.org/10.1111/febs.15219](https://doi.org/10.1111/febs.15219)





# ChIP-Seq Results Folders/Files and Report




# Folder structure and file types for QC and peak calling results returned by CAB

- Four folders: Coverage, Peaks, Stats, SRM#\_LabIDs
- Two files: HTML report, PDF report
- Within SRM#\_LabIDs: Coverage, Peaks, Stats
- SRM#\_LabIDs contains files labeled by ChIP targets
- Coverage folders will contain BigWig files
  - Use these to visualize data in IGV, etc.
- Peaks folders will contain peak calling results
  - Use these to visualize peaks in IGV, etc.
  - Contain “filter” version of each peak files
  - Filter files remove ENCODE “blacklist” regions
- Stats folders will contain four QC-related files
  - Mapping metrics, QC metrics
  - Cross-correlation analysis table
  - Cross-correlation plot



# Common options for visualizing data

- [Integrative Genomics Viewer](#) (IGV) - download desktop app or use web browser
- [UCSC Genome Browser](#) – select reference genome, upload your files (custom tracks)
- [St. Jude ProteinPaint](#) – select reference genome, upload your files
- Files to view often includes BigWig, broadPeak, narrowPeak, bam files with expression data (e.g. RNA-Seq)



Adapted from [hbctraining.github.io/Intro-to-ChIPseq](https://hbctraining.github.io/Intro-to-ChIPseq)



# Example ChIP-Seq report from CAB

---


- Interactive HTML report, static PDF report
- HTML, PDF examples of reports available on CAB's wiki page
- Walk through HTML example
- Download HTML to follow along [here](#)



# **Additional ChIP-Seq Analyses from CAB**

# Examples of other analyses that can be requested from CAB: Differential binding

- Additional analyses are considered customized or collaborative
- Require additional SRM requests
- For example:
  - Differential binding site analysis
  - Chromatin state assignment
- Differential binding site analysis can identify differences in enrichment peaks across conditions



Select the data type\*

ChIPseq/Cut-and-Run

Select a bioinformatics analysis from the list below. If it's not on our list, please choose "Other" and describe the analysis in the box to the side.

Bioinformatics analysis\*

- ☐ A1a. QC and Peak calling
- ☐ A1b. Peak annotation (ChIPseq)
- ☐ A2. Differential binding analysis (ChIP)
- ☐ A3. Chromatin state assignment (chromHMM)
- ☐ A14. Competitive Mapping (for spike-in or de-contamination)
- ☐ A15a. GO enrichment analysis (ChIP WGBS)
- ☐ A16a. DNA motif analysis (ChIPseq/WGBS)
- ☐ A17. Super-enhancer analysis and Circular regulatory circuitry (H3K27ac)
- ☐ A18. Co-localization analysis
- ☐ A26. Spike-in or SpikeInFree normalization
- ☐ U2. Mapping (ChIPseq)
- ☐ Other



done

Adapted from [doi.org/10.3389/fgene.2015.00169](https://doi.org/10.3389/fgene.2015.00169)



# Examples of other analyses that can be requested from CAB: Chromatin state assignment

- Chromatin state assignment uses ChromHMM to predict state
  - Based on Hidden Markov Model (HMM), models presence or absence of chromatin marks, annotates genome
- Hidden Markov Models are probabilistic, predict “outcomes” based on observable parameters
- Includes “hidden” states that influence outcome but aren’t “observable”
- In this case, “outcome” is chromatin state, predicted based on observed chromatin mark patterns



Credit: [epicypher.com](http://epicypher.com)





**Questions?**